



PRODUCTS AND EXPERTISE FOR DISTRIBUTORS AND HEALTHCARE OPERATORS

Positioning Systems for the Operating Room

www.biomatrix.it





Summary

Introduction	4
Vacuum mattresses Bi-VACUUM	8
Bi-Pad disposable antislip mat for Trendelenburg operations	12
Prone facial positioner	14
Air positioning pillow	15
Prone positioners	16
Octo-Pos Positioners	17
Gel positioners	18
Disposable positioners	28
Multi-purpose polyurethane positioners	32
Hygiene	36

Operational uses

○ Supine position



It gives access to almost all surgical sites. It requires the head to be aligned with the rest of the body.

ATTENTION! Do not let feet protrude from the table to avoid compression or injuries of the Peroneal nerve or the Achilles tendon.

○ Gynecological or Lithotom Position



This position requires the intervention of two people who raise the lower limbs simultaneously. If the patient is awake the medical staff could ask for collaboration, obtaining the desired position gradually. The same care should be taken in returning to the primary position.

ATTENTION! The perineal area must protrude from the tabletop by about 3–5 cm, in order to carry out the intervention.

○ Lateral Decubitus



Position used for thoracic surgery and for kidney surgery. Being a difficult position to obtain, the use of a comfortable support surface is essential.

ATTENTION! The sling (lower arm) must not exceed a 90° angle to the body. The flexion is applied differently depending on the intervention to be performed.

○ Prone position



This position allows surgical access to the spine and the back of the head.

ATTENTION! The head must be aligned with the rest of the body and the face positioned on a suitable device that allows the intubated area to be left free.

○ Recommended positioners: **Supine position**

Head



Upper Limbs




Sacral area



Thighs and Knees



Legs and Feet



○ Recommended positioners: **Gynecological or Lithotom Position**


Head



Upper Limbs




Shoulder area



Loin-sacral area



Lower Limbs



○ Recommended positioners: **Lateral Decubitus**

Upper Limbs



Head



Lower Limbs



Pelvis



Knees



Feet and Ankles



○ Recommended positioners: **Prone position**


Head




Upper Limbs



Thorax



Pelvis




Lower Limbs




Positioners for the operating Room

POSITIONING SYSTEMS FOR THE OPERATING ROOM




○ Type of intervention or Surgery technique

After studying the case, the surgeon, and the team, define the optimal position considering the specific patient's conditions.




○ Anatomical area to be treated

Depending on the localization of the intervention, it is possible to prefer various types of positioners. Size, conformation, and type of positioner can be determined not only by the district to be treated, but also by the anthropometric characteristics of the patient. It is also essential to consider the potential pressure damage to areas not directly affected by the surgical activity.



○ Type of material and positioner

Finally, the duration of the intervention, the level of comfort required, the type of practices and the tools used to verify the compatibility of materials are considered.



Bi-VACUUM

Vacuum Mattresses

Ideal for long-lasting operations and for interventions that require the exposure of anatomical areas that cannot be reached in conventional postures. Particularly recommended for elderly patients or people with compromised physical and motor conditions.

Depression can be carried out even during surgery

The process of introducing and vacuuming the air inside the mattress can also be carried out during the intervention, thanks to the non-return valve. If deemed necessary by the surgeon, the sedated patient can be repositioned without Bi-VACUUM being removed indeed. The patient can also be subjected to radiological examinations thanks to Bi-VACUUM radiolucency.

Use it as a transfer device


Once the intervention is completed, the mattress can be positioned as if it were a transfer sheet and depressed to give it the necessary rigidity to use it as an aid to transfer the patient.

Outer layer

The outer part, composed of a mix of PVC and Polyurethane, is characterized by a velvety effect in contact with the patient, which gives comfort and prevents the onset of any skin irritation or lesion.


Inner composition

Its internal composition, consisting of expanded polystyrene microspheres, stabilizes the patient optimally, adapting ergonomically to his anatomy and position taken.





POSITIONING SYSTEMS FOR THE OPERATING ROOM



CPC valve
Non-return valve.



● **Bcode**

Find on-line info on this category right now just by shooting this code

Bi-VACUUM

Bi-VACUUM brings numerous advantages to the work of operators, ensuring the highest standards of comfort and safety for patients:

Greater stability

Bi-VACUUM vacuum positioners conform perfectly and firmly to the patient's body and guarantee a safer positioning than that obtained with other systems which are simply placed against the patient.

Access to the limbs

The Bi-VACUUM vacuum positioners allow you to immobilize the patient in positions that are difficult to reach with other systems. The Bi-VACUUM vacuum positioners also allow to maneuver the patient's legs and arms without having to reposition.

Excellent access to the operating field


Thanks to the Bi-VACUUM positioners, the patient can be positioned exactly how the surgeon wishes for utmost preparation of the operating field.

Time savings


Using the Bi-VACUUM vacuum positioners, the patient is easily placed on the operating table and, if necessary, it can be quickly repositioned during the intervention saving valuable time to the entire staff.

Standard Mattresses


+ NEW




SO058006




SO058003




SO058004



SO058002



SO058005




SO058001

Code	Size
SO058006	100x200 cm
SO058003	75x200 cm
SO058004	65x200 cm


Code	Size
SO058002	75x125 cm
SO058005	100x100 cm
SO058001	75x100 cm

Repair Kit AVAILABLE
SO05W00 repair kit for Bi-Vacuum depression mat


Shaped mattresses




SO059004



SO059003



SO059002




SO059001


Code	Size
SO059004	75/100x125 cm
SO059003	80/95x105 CM

Code	Size
SO059002	75x100 cm
SO059001	35x40 cm


Standard mattresses with 4 chambers/ Shaped mattresses with 4 chambers




SO05900S3



SO05900S2



SO05900S1



SO05800S1

Code	Size
SO05900S3	75/100x125 cm
SO05900S2	100x100 cm

Code	Size
SO05900S1	80x100 cm
SO05800S1	80x200 cm

Bi-Pad

Disposable antislip mat for Trendelenburg operations

It fits perfectly to the patient's body, avoiding slipping and favoring its anchoring to the operating table during interventions that require the position of Trendelenburg and reverse Trendelenburg (gynecological surgery, colorectal, urological); performed both laparoscopically and robotics.

Made of slow memory viscoelastic polyurethane, this material performs an important anti-bedsores function, protecting the patient from skin pressure injuries.

Its disposable feature also ensures minimizing the risks of cross contamination.

Bi-Pad

+ NEW

STANDARD Dimensions

Length: 90 cm

Width: 52 cm

Total height: 3 cm



Non-slip
base

SO02815DVG52

Bi-Pad disposable pad

SO02815S

Bi-Pad with silicone

SO02815V

Bi-Pad with straps and spunlace sheet

SO02815V.I

Bi-Pad with cut-out and spunlace sheet

SO02815V.75

Bi-Pad with straps and spunlace sheet
cm 75x50x3

SO02815VKIT

Bi-Pad with straps, spunlace sheet
and thorax belts

SO02815VKIT.I

Bi-Pad with cut-out, straps, spunlace sheet
and thorax belts


SO02815VL

Bi-Pad long with straps, spunlace sheet, cm 180x52x3

SO02815VKITL

Bi-Pad long with straps, spunlace sheet and thorax
belts, cm 180x52x3

Chance to create customized KIT



Example of use of the non-slip mattress in the operating room (in combination with SO02800 - Disposable donut headrest; SO059004 - Shaped vacuum mattress).


Spunlace tubular sheet

The tubular Spunlace sheet can be used to restrain the patient's arms, ensuring greater stability on the operating table during interventions in Trendelenburg.

SO02W815S

+ NEW

Spunlace tubular sheet



Technical features

Material: Spunlace fabric

Dimension: 220x30 cm

Packaging: 10 pieces



Example of use of the Spunlace sheet (SO02W815S) in combination with the non-slip mat (SO02815V).

Disposable Prone Facial Positioner and multi-purpose Mirror

Multipurpose mirror for disposable prone facial positioner.



SO02816 - SO02W816

+ NEW

Ideal where there is a need to keep the patient in prone position.

The high-density polyurethane foam that composes it, gives comfort to the patient, avoiding unpleasant discomfort and injury during and after prolonged prone positioning of the face. The wide support surface on the operating table, gives stability to the head during interventions, while ensuring comfort and safety.

Thanks to the application of a multi-purpose acrylic mirror, which is placed at the base of the face positioner, the surgeon and the medical team can constantly monitor the patient's condition during surgery, whether the patient is completely or partially sedated. The synthetic material of the mirror makes it extremely light and at the same time safe (SO02W816).

The face positioner is also used in Intensive Care Unit for patients with respiratory dysfunctions who have the need to assume the prone position to improve the lung capacity: in this situation and in all cases in which the patient is prone, conscious and not sedated, the positioner does not need the mirror.



Example of use of the Prone Facial Positioner in the operating room, with mirror for patient monitoring (in combination with SO03800 - Prone positioner; SO02813 - Disposable shaped mat and CD01608 - Multipurpose cylindrical cushion in polyurethane).

Air positioning pillow

Inflatable positioning cushion designed to support the patient's shoulder during surgical procedures and examinations.



SO04800

+ NEW

Air positioning pillow

The device can be adjusted manually, depending on the necessity, ensuring good weight distribution and a noticeable improvement in patient comfort.

Its use allows to lift and turn patients in a quick, simple, ergonomic, hygienic, and safe way.

The device is useful during total hip replacement, thoracotomy, shoulder arthroscopy and kidney surgery.

Its use is recommended with slightly sedated patients undergoing epidural anesthesia and patients undergoing long-term interventions.

The air positioning cushion can be easily cleaned with a cloth dipped in any mild detergent or alcohol 70% isopropyl or 0.5% bleach solution or with a disinfectant wipe. It is recommended to make sure the device is completely dry before using it again.



SO04800 56X25 cm

SO04801 63X25 cm

SO04802 50X25 cm

SO04803 30X20 cm

Technical features

Material: entirely of nylon-polyurethane.

Supplied with inflation tube.

Permeable to X-rays. Latex free.



On request, it is possible to substitute inflating valve with CPC valve

Prone positioners

They offer a versatile and efficient solution for the support and positioning of the patient, in prone position, during medical procedures.

SO03800 - SO03800S

+ NEW

Prone positioner

Its structure allows the uniform distribution of the pressure even with patients of very high body weight. Reduces pressure in the abdominal area, and while used during interventions, reduces abdominal hypertension bringing some advantages, like better ventilation and bleeding reduction.

The prone positioner is made of high-density fire-retardant polyurethane foam and covered with a removable cover in fireproof Skay.

Technical features

The inner hole has a trapezoidal shape, in the area in contact with the patient. This particular feature makes it more comfortable and safe.

SO03800 80X50X12 cm

SO03800S 80X46X10 cm

COVERING IN
FIREPROOF
SKAY



Chest positioning cushion

Made internally in flame retardant polyurethane viscoelastic, density VE50 / FR-VE850385 and covered in fireproof eco-leather. The materials used are all Italian.

It is recommended to use the covered device with a disposable TNT towel to ensure one greater hygiene.

Item can be realized in different dimensions.

Dimensions

Length: 30 cm

Depth: 50 cm

Height: 17,10 cm

SO03804

+ NEW





SO03800N

+ NEW

Octo-Pos

Increased comfort thanks to the new rounded design

Thanks to the new design, the positioner offers rounded shapes and corners, especially on the critical contact points. This way, pressure damage is avoided, particularly during long interventions.

The inside has a 45° beveling

Total flexibility thanks to the modular additional spacers

As standard, the positioner is equipped with removable modules of 5cm thickness which allow the device to be optimally adjusted to the patients. For an even more precise sizing and in order to meet any need of the operating room, 2.5cm spacer kits are available.

Made of foam coated with polyurethane water repellent, single component, breathable, high strength, formaldehyde-free, bio-rigidly degradable resin; it contains no toxic additives.

SO03800N 76X45X11,5 cm






Clavicles and thorax
beveling



One single device
also for oversized
patients

Standard
spacers

Optional
spacers



Gel Positioners for the operating Room

The person undergoing surgical surgery, partially or totally sedated, may be forced in the same position for hours, sufficient time for the onset of pressure injuries.

The Gel Positioners are compact and thanks to their physical characteristics and their density, they can adapt to the body surface, allowing a comfortable support of the patient on the operating table.

Among the properties we recognize the following:

- water repellent capacity,
- their hardness does not change over time,
- if damaged by cuts gel does not come out,
- easily sanitizable,
- isothermal,
- can be used in combination with warming or cooling mattresses.

Wanting to increase comfort and prevention of skin lesions the gel pads can be combined with the foam positioners.



● **Bcode**

Find on-line info on this category right now just by shooting this code ice

Line of Gel Positioners



POSITIONING SYSTEMS FOR THE OPERATING ROOM



SO01800

Closed ring for head

- SO01800S - 9x4x2 cm
- SO01800M - 14x6x3,5 cm
- SO01800L - 20x7,5x4,5 cm
- SO01800XL - 20x7,5x7,5 cm



SO01801

Open ring for head

- SO01801S - 9x4x2 cm
- SO01801M - 14x6x3,5 cm
- SO01801L - 20x7,5x4,5 cm
- SO01801XL - 20x7,5x7,5 cm



SO01802

Horseshoe ring for head

- SO01802S - 9x4x2 cm
- SO01802M - 14x6x3,5 cm
- SO01802L - 20x7,5x4,5 cm

Gel Positioners

for the Operating Room

POSITIONING SYSTEMS FOR THE OPERATING ROOM



SO01804

Prone headrest


SO01804L
28x23x14,5 cm
SO01804S
23x19x13 cm



SO01803

Face pad

SO01803
16x6x1,6 cm



SO01805

Contoured head rest

SO01805
19x19x6 cm



SO01806

Head pad

SO01806
27x23x7 cm



SO01807

Supine head rest


SO01807
27x24x9 cm



SO01849

Prone headrest 4 ways

SO01849
23x20x11 cm



SO01849-1

Prone headrest 2 ways

SO01849-1
23x20x6 cm

Gel Positioners

for the Operating Room



SO01808

Armboard pad


- SO01808S
40x12,5x2 cm
- SO01808L
60x12,5x2 cm



SO01809

Anatomic armboard pad

- SO01809S - 40x12,5x6 cm
- SO01809M - 50x12,5x6 cm
- SO01809L - 60x12,5x6 cm
- SO01809XL - 63x12,5x6 cm
- SO01810 - 52x15x5,5 cm



SO01824

Operating table pad with foam


- SO01824A - 16x8x6 cm
- SO01824B - 70x20x8 cm
- SO01824C - 30x15x10 cm
- SO01824D - 50x15x10 cm
- SO01824E - 60x15x10 cm



SO01811

Ulnar nerve protector

- SO01811S
38x12x0,8 cm
- SO01811L
46x21x0,8 cm



SO01820

Lateral positioner

- SO01820
72x50x16 cm



SO01812

Fixing belts

- SO01812S
33x7x1,5 cm
- SO01812L
36x7x1,5 cm

Gel Positioners

for the Operating Room



SO01814

Shaped chest roll


- SO01814S
30x16x13,9 cm
- SO01814M
40x16x13,9 cm
- SO01814L
50x16x13,9 cm



SO01815

Chest positioning set


- SO01815SX
71/58x20x7,6 cm
- SO01815DX
71/58x20x7,6 cm



SO01815

Thyroid Gland positioner for surgical interventions

- SO01851
51x40x14 cm



SO01825

Shoulder pad with Velcro

- SO01825
5x13x2,5 cm



SO01813

Chest roll – Positioning roll with flat base

Dimensions • Adults

Dimensions • Children

- | | |
|-------------------------|--------------------------|
| SO01813A - 30x15x10 cm | SO01813H - 15x5x5 cm |
| SO01813B - 40x15x10 cm | SO01813I - 20x5x5 cm |
| SO01813C - 50x15x10 cm | SO01813J - 25x5x5 cm |
| SO01813D - 30x10x8,5 cm | SO01813K - 30x5x5 cm |
| SO01813E - 40x10x8,5 cm | SO01813L - 15x7,5x7,5 cm |
| SO01813F - 50x10x8,5 cm | SO01813M - 20x7,5x7,5 cm |
| SO01813G - 50x18x15 cm | SO01813N - 25x7,5x7,5 cm |
| | SO01813O - 30x7,5x7,5 cm |

Gel Positioners

for the Operating Room



SO01828

Dome positioner

SO01828
35x15x7,5 cm



SO01816

Hip positioner

SO01816
14x20x1,3 cm



SO01817-18

Hip pad with hole for pin


SO01817
18x18x1,3 cm
SO01818
33x30x0,8 cm



SO01819

Universal patient positioner

SO01819A - 16x8x4,5 cm
SO01819B - 33x7x5 cm
SO01819C - 46x21x0,8 cm
SO01819D - 36x12,5x7 cm
SO01819E - 25x12,5x5 cm
SO01819F - 38x12,5x6,3 cm
SO01819G - 50x12,5x7 cm



SO01821

Relton frame set

SO01821
19x17,5x10,5 cm

Gel Positioners for the operating Room

Wilson Frame

SO01822C - SO01823

Set telaio di Wilson

It is used to keep patients in a flexed position during lumbar and spinal procedures (lumbar laminectomy, decompressions, spinal surgery).

The gel pads protect the patient and grant a more comfortable position. The protections guarantee a low pressure on the patient skin.

SO01823

Wilson frame 60x49x15,5 cm

SO01822C

Lateral long pad. 72,5x17,5x2 cm

Non flared pad 18,5x12x2 cm

Flared pad 18,5x12x2 cm



Wilson frame pads

SO01822

They are fixed to the Wilson frame thanks to appropriate belts; the gel coating grants a more comfortable position.


The protections guarantee a low pressure on the patient skin.

SO01822A

75x25x2,5 cm

SO01822B

75x25x1 cm





SO01826A

Heel supports

SO01826A
20x11x7 cm



SO01826B

Heel protectors pads

SO01826B
17x11x5 cm



SO01827A

Heel pads with straps

SO01827A
17x18x1,3 cm



SO01827B

Heel pads with Velcro


SO01827B
17x23x1,3 cm



SO01832

Boot stirrup pad

SO01832
71x41/30x1 cm



SO01833

Traction boot pad

SO01833
37x24x0,5 cm

Operating table pad

Operating table pads protect the whole operating table surface and they are recommended for every position.

These pads are soft, comfortable and versatile. They distribute the patient's weight on the entire surface.


SO01837 - SO01839



- SO01837A - 23x50x1 cm
- SO01837B - 99x50x1 cm
- SO01837C - 50x50x1 cm
- SO01837C-1 - 50x50x1,5 cm
- SO01837D - 50x25x1 cm
- SO01837D-1 - 50x25x1,5 cm
- SO01837E - 10x38x1 cm
- SO01837F - 42x15x1 cm
- SO01837G - 10x61x1
- SO01837H - 50x45x1 cm
- SO01837H-1 - 50x45x1,5 cm

- SO01837I - 193x86x1 cm
- SO01837J - 183x50x1 cm
- SO01837J-1 - 183x50x1,5 cm
- SO01837K - 117x50x1 cm
- SO01837K-1 - 117x50x1,5 cm
- SO01837L - 160x43x1 cm
- SO01837M - 15x35x1 cm
- SO01837M-1 - 15x35x1,5 cm
- SO01837N - 70x50x1 cm
- SO01837N-1 - 70x50x1,5 cm
- SO01837O-1 - 70x45x1,5 cm

- SO01837P-1 - 15x50x1,5 cm
- SO01839A - 46x33x0,5 cm
- SO01839A-1 - 46x33x1 cm
- SO01839B - 22x17x1,2 cm
- SO01839C - 60x60x2,5 cm
- SO01839D - 90x50x2,5 cm
- SO01839E - 75x25x2 cm
- SO01839F - 10x15x2 cm
- SO01839G - 15x66x1,3 cm
- SO01839H - 70x50x1,6 cm
- SO01839I - 100x50x1,6 cm



SO01836

Embossed operating table pad


- SO01836A - 50x50x1,6 cm
- SO01836B - 115x50x1,6 cm
- SO01836C - 40x23x1,6 cm
- SO01836D - 91x45x1,6 cm



SO01838

Operating table pad with notch

- SO01838A - 117x50x1 cm
- SO01838A-1 - 117x50x1,5 cm
- SO01838B - 100x50x1 cm
- SO01838B-1 - 100x50x1,5 cm
- SO01838C - 50x50x1 cm
- SO01838C-1 - 50x50x1,5 cm



SO01847

Sacral pad

- SO01847A cm 40x40x1,6
- SO01847B cm 46x40x1,6
- SO01847C cm 46x46x1,6
- SO01847D cm 50x46x1,6
- SO01847E cm 50x50x1,6

Gel Positioners

for the Operating Room



SO01831


Lloyd support protection

SO01831A - 31x23x1 cm

SO01831B - 40x25x1 cm

SO01831C - 52x30x1 cm

SO01831D - 55x52x1 cm




SO01829

Graduated positioner

SO01829

23x36x3,8 cm



SO01830

Cot positioner

SO01830

30,5x66x1,27 cm



SO01834

Lithotomy cover

SO01834

71x12x0,8 cm

(SET of 2 piece)



SO01835

Ankle stirrups

SO01835A

31x4x1 cm

(SET of 2 piece)

SO01835B

44x4,5x1 cm

(SET of 2 piece)



SO01840

Universal patient positioner

SO01840

50,5x8,5x0,8 cm

Disposable positioners for the operating Room

The disposable Bi-FOAM positioners for the operating room are designed to guarantee patient utmost comfort e prevent the risk of contamination during surgical procedures.



Bcode

Find on-line info on this category right now just by shooting this code

Polyurethane devices are a worth alternative to the reusable positioners, sometimes they can be combined with gel pads to increase patient comfort on the operating table above all during long-term interventions

They optimize the final positioning indeed, which would not be obtained with gel devices only.

The Bi-FOAM also create excellent non-slip grip in positions that require maximum inclination. Recommended in dirty interventions in order not to damage or pollute the gel pads, which, due to accidental cuts or damages could be infected, or polluted by material that is difficult to expel.

The disposable positioners are made of polyurethane foam, elastic and non-deformable, 25kgm³ +/- 5%, non-toxic and latex free.

Code	Size	Dimintions	Shape
SO02800	U	ø 23 x 5 cm	
SO02803	U	23x23x6/10	
SO02802	U	23x22x5/10	
SO02805	U	13,5x61x7v	
SO02815B8	U	90x8 cm	
SO02807	U	42x26,5x14	
SO02808	U	17x18x3,5 cm	
SO02814	U	7,6x152,4 cm	
SO02806	U	38x15,5x3,5 cm	
SO02813	U	183x51x5 cm	



SO02800

Donut headrest



SO02802

Adult head positioner with hole



SO02803


Slotted headrest with hole

SO02800 diametro cm 23x5

SO02802 23x22x5/10

SO02803 23x23x6/10 cm

Disposable positioners for the Operating Room



SO02806

Shaped pads

SO02806 38x15,5x3,5 cm



SO02805

Standard laminectomy arm cradle


SO02805
13,5x61x7 cm



SO02807

Anatomic leg positioner

SO02807
42x26,5x14 cm



SO02808

Shaped heel protector

SO02808
17x18x3,5 cm



SO02811

Arm Support

SO02811
55x35x2,5 cm



SO02813

Shaped pad

It offers a low-cost and effective solution for patient support during all types of surgeries. Can be used in combination with other positioners also and especially during operations where the patient is semi-conscious.

SO02813
183x51x5 cm

Disposable positioners for the Operating Room



SO02820A - SO02820B

+ NEW

Disposable Cylindrical chest roll

Disposable cylindrical chest roll, made of elastic and non-deformable polyurethane foam 25 kgm³, non-toxic and latex free, non-fireproof.

SO02820A

Ø 15cm - length cm 40

SO02820B


Ø 15cm - length cm 60

Technical data

Material:

T25 MKS polyurethane yellow color (Italian production)

Packaging: Individually wrapped and vacuum packed



SO02823

+ NEW

Chest half roll positioner in PU

Disposable positioner with a semi-cylindrical roller, made of elastic and non-deformable expanded polyurethane with a surface in memory foam, non-toxic, non-fireproof.

SO02823

20x23x9.5 cm (± 2%)

Technical data

Material:

Non-Fireproof polyurethane foam. The fireproof model is available upon request. Disposable - Non sterile

Packaging: Individually packed in a transparent plastic bag.



SO02830

+ NEW

Adjustable Lateral Head Positioner

Disposable PU positioner to accommodate patients' head in case of surgery requiring lateral positioning. Made of PU, the upper top layer is in viscoelastic memory foam to counteract any effect of pressure on the patient's face and avoid decubitus damage. The structure of the positioner is composed of easily adjustable layers that allow different height adjustment options to obtain correct spinal alignment and a central hole for the patient's ear. The ergonomic design also allows use for prolonged lateral interventions.

SO02830

28x20x17/25 cm

Technical data


Materiale:

Non-Fireproof polyurethane foam. The fireproof model is available upon request. Disposable - Non sterile

Packaging: individually wrapped in transparent plastic bag

Disposable positioners

for the Operating Room



SO02814

Limb holder with Velcro

It offers an low-cost and efficient solution for light patient restraint during medical procedures. Limb support made of polyurethane foam e polyamide, thickness mm 6 +/- 0.3, with Velcro closure.

SO02814
64X5 cm (20 pcs/box)

Technical data

Material:

Polyurethane foam 54%
Polyamide 46%



SO02815B8

+ NEW

Fixing thorax belts

Made of polyurethane foam and polyamide, thickness mm 6 +/- 0.3, with Velcro closure.

SO02815B8
140X8 cm (LxH)



SO02819A - SO02819B

+ NEW

Velcro bands for legs

SO02819A
Velcro Bands for legs L-size
cm 220x10 (64 pz/box)

SO02819B
Velcro Bands for legs M-size
cm 150x10 (80 pz/box)



Multi-purpose polyurethane positioners for the operating room and post-operative hospitalization

Line of positioners made of polyurethane, with a load bearing structure, coated with a non-toxic waterproof paint. Like the Bi-FOAM range in lightness, breathability and in ease of use, they differ in the particular film that covers them allowing an easy sanitization and making them reusable.

They are designed to adapt and support the different parts of the patient's body, both during surgery and in the immediately following phase, i.e. during post-operative recovery.

Cleaning can be carried out with disinfectants diluted up to 5% even alcohol-based, or with a cloth soaked in water (up to 60 ° C) and neutral detergent.



CD01600

Donut cushion

It protects and supports the head, ears, the neck, and face during operations.

CD01600S

ø 20 cm

CD01600L

ø 24 cm



CD01606

Retractor

Used to ensure an assisted and comfortable opening of the legs.

CD01606S

15x9x7,5/13 cm

CD01606L

15x16x14/20 cm




CD01610

Heel protector for clubfoot

Used in the presence of clubfoot.

CD01610

23x46 cm




CD01601

Headrest

It protects and supports the head during surgical procedures.

CD01601
27x24x9 cm



CD01602

Trapezoidal positioning cushion

Suitable for tonsils operations, intubation, and ventral positions.

CD01602S
30x25x10/5 cm
CD01602M
40x25x6/2,5 cm
CD01602L
40x25x10/5 cm



CD01603

Triangular positioning cushion

Suitable for dorsal anesthesia.

CD01603S
30x20x10 cm
CD01603L
40x25x18 cm



CD01604

Triangular positioner

Suitable for lateral decubitus up to 30°.

CD01604S
55x25x18 cm
CD01604L
70x25x18 cm



CD01608

Cylindrical cushion

Suitable for the protection of knees and pelvis.

CD01608S
ø 40x15 cm
CD01608L
ø 60x15 cm



CD01609

Armrest with handle

Used to position and relieve arm and shoulder.

CD01609
50x12x15 cm

Polyurethane positioners


for the operating room

Positioners made of load-bearing polyurethane or cold-molded polyurethane, free of protective film, but supplied with a removable cover.


Depending on the trauma and the type of surgery reported by the patient, aids with covers will be recommended: in cotton, or in breathable bielastic polyurethane.



CD01605D
40x47x20/14 cm
Belts:8 cm

 Rehabilitation

 Intensive Care Unit

 Long-term care

Disposable inguinal retractor

Retractor in polyurethane foam with lateral grooves, available in different densities depending on the destination of use: single use, if used in the operating room, or single patient, if used during the recovery from hip surgery or rehabilitation from other trauma.

Accessories

It can be supplied with one or two adjustable fixing straps with Velcro closure, to immobilize the patient's lower limbs and / or to fix the retractor to the operating table.

 CD01605D

+ NEW



 CD01605

Inguinal retractor with bielastic polyurethane cover

Shaped retractor, favours low pressure in the support areas.

CD01605
40x47x20/14



 CD01607

Inguinal retractor with cotton cover

Ideal positioning device to support the lower limbs in the postoperative period.

CD01607S
44x27x13 cm
CD01607L
50x35x20 cm



 CD01610PO

Post-surgery splint

Used in orthopedics and traumatology to restraint and immobilize the lower limb after fracture or surgery.

CD01610PO
73x35x53 cm



Hygiene



Hygiene

A selection of reliable, safe Hygiene Devices, produced with high quality and technical materials.

Due to the inflow of patients, technical health personnel and frequent contact with external agents, hospital departments can act as a vehicle for the proliferation of bacterial and viral contamination and infections.

The hygiene of medical staff and patients is therefore decisive to break down as much as possible the risk of contamination that could compromise the health of the inpatients.

● Bcode

Find on-line info on this category right now just by shooting this code



○ Oral hygiene of the patient

La cavità orale è fonte di infezioni quali polmoniti nosocomiali e da ventilazione associata (VAP).

The presence of dental plaque, mucus and secretions facilitates bacterial colonization of the area and the proliferation of pathogens, increasing the patient's stay in the intensive care unit.

To prevent the onset of such pathologies, it is appropriate to develop and implement a complete hygiene protocol that may include the use of disinfectants, the use of appropriate safeguards and hydration of the oral cavity, to protect the patient's health in the best possible way, while also allowing the hospital to save money.

○ Hygiene of the hands for operators and inpatients

It is fundamental to grant Hands' Hygiene of both healthcare personnel and patients during the preoperative procedures, to reach excellent hygienic conditions and prevent any risk of infection.



Oral Swab

Ideal in:



Intensive Care Unit

Three different types of oral hygiene swabs for patients, partially collaborating or intubated.

The operator must guarantee oral hygiene assistance to the patient if the latter is conscious or partially conscious; in more serious cases, the operator must personally perform all operations.

Cleaning Swab



IG011000-IG011000B

Cleaning Swab in addition to removing residues present in the oral cavity stimulates the gingival tissue. Equipped with a soft foam head with marked ridges, it also allows interdental cleaning eliminating the part of devitalized tissue. It can also be used for the application of soothing creams.

Available with (IG011000B) and without bicarbonate (IG011000).

Aspiration Swab



IG011001-IG011001B

Aspiration Swab allows removing debris and secretions during the different oral cavity cleaning operations. Thanks to the soft foam head it cleans gently, stimulating the teeth and the gum tissue. It is equipped with an ergonomic grip for the thumb, ensuring absolute control of the device with the use of a single hand. The end of the attachment, being conical, can adapt to any type of vacuum cleaner.

Available with (IG011001B) and without bicarbonate (IG011001).

Brush Swab



IG011002-IG011002B

Brush Swab allows you to remove secretions and rests in the oral cavity during cleaning procedures. Thanks to soft foam head gently cleans teeth stimulating the gum tissue, moreover the presence of bristles effectively removes dental plaque. It is equipped with an ergonomic thumb grip, ensuring control of the device with the use of only one hand. The end of the attachment, being conical, can adapt to any type of vacuum cleaner.

Available with (IG011001B) and without bicarbonate (IG011002).

Technical features

Box: 1000 pcs

Packaging: singolarmente

Use: One use only

Material: Soft foam head with marked ridges

Technical features

Box: 800 pcs

Packaging: singolarmente

Use: One use only

Material: Latex free

It is equipped with an ergonomic grip for the thumb

Technical features

Box: 600 pcs

Packaging: singolarmente

Use: One use only

Material: Latex free

Soft Foam head Bristles



HYGIENE

Hand hygiene for operators and patients

Surgical preparation of the hands represents a critical moment for the safety of health care.

It aims to reduce the release of cutaneous bacteria from the hands of the surgical team for the whole duration of the procedure, so that in case of hypothetical undetected breakage of the surgical glove, no bacteria contaminate the operating field.

Ideal in:



Operating room



IG011010

Surgical scrub brush

Designed for comfort and adaptability. Gentle on the skin, avoids irritation, reducing the risk of infections. Soft bristles for nails and knuckles and sponge for hands and forearms skin.

Equipped with a nail cleaner. Sterile product.

Technical features

Box: 450 pcs (15 dispenser of 30 units)

Dispenser: 100 pcs

Use: One use only

Available in versions:

without disinfectants (**IG011010**)

with chlorhexidine (**IG011010C**)





POSITIONING SYSTEMS FOR THE OPERATING ROOM

BioMatrix

Via Maffini, 23 - 29010 San Giuliano di Castelvetro (PC)
Tel. 0523.826619 Fax. 0523826728
e-mail: info@biomatrix.it

All right reserved. Changes may be made without prior notice.
BioMatrix S.r.l has not to be considered responsible for damages
resulting from the lack or inaccuracy of the information reported here.

www.biomatrix.it

